

A Level Chemistry B (Salters)
H433/03 Practical skills in chemistry
Sample Practical Insert

Date – Morning/Afternoon

Time allowed: 1 hour 30 minutes



INFORMATION FOR CANDIDATES

- This document consists of 4 pages. Any blank pages are indicated.

INSTRUCTIONS TO EXAMS OFFICER/INVIGILATOR

- Do not send this insert for marking; it should be retained in the centre or destroyed.

Iron and manganese in paper clips

A student describes below a project to find the amount of iron and manganese in some paper clips:

- To find the amount of iron in the paper clips, I decided to use a titration with potassium manganate(VII) solution.
- To find the amount of manganese in the paper clips, I found out that I could oxidise it in solution to potassium manganate(VII) and then use a colorimeter because the intensity of the purple colour I get depends on its concentration of MnO_4^- ions.

Part 1: Determination of iron in paper clips

I found the following method in an old book and decided to use it to find the amount of iron in the paper clips.


Determination of the percentage of iron in iron wire

Weigh out accurately about 1.4 g of iron wire and transfer it to a conical flask containing 25 cm³ of dilute sulfuric acid and a few cm³ of concentrated sulfuric acid to accelerate the reaction. Fit the flask with a rubber bung containing a short length of capillary tubing.


Warm the flask carefully to maintain a steady reaction and, when all the iron has reacted leaving only particles of carbon, cool the flask. Transfer the solution quantitatively to a 250 cm³ volumetric flask and make up to the mark with dilute sulfuric acid and water. Pipette 25 cm³ of this solution into a conical flask, add about 25 cm³ of dilute sulfuric acid and titrate with standard potassium manganate(VII) solution (about 0.02 mol dm⁻³).

[Reference: STARK, J G (1971): *Titrimetric analysis for A & S levels SI Edition* London, John Murray, 27]

The reaction of iron with sulfuric acid is:



The half-equations for the titration reaction are:



I had to ask what a 'capillary' tube was and was told that it is a glass tube with a small internal diameter.

These are my results:

Mass of paper clips added to the conical flask = 1.28 g

Concentration of manganate(VII) solution = 0.0200 mol dm⁻³

	Titration 1	Titration 2	Titration 3	Titration 4	Titration 5
Final burette reading / cm ³	22.90	45.40	22.55	43.05	24.00
Initial burette reading / cm ³	0.00	22.90	0.00	20.95	1.55

Part 2: Determination of manganese in paper clips

I found a website that said that steel normally contains between 0.1 and 0.4% manganese. I was given the following worksheet that I could follow to find the amount of manganese in the paper clips. It involves reacting pieces of paper clip with nitric acid to produce a solution containing Mn^{2+} ions. The Mn^{2+} are then oxidised with potassium iodate(VII) to MnO_4^- ions.

Method for the determination of manganese in paper clips:

- 1 Weigh accurately about 0.25 g of cut-up paper clip.
- 2 Put it into approximately 70 cm³ of 2.0 mol dm⁻³ nitric(V) acid in a beaker.
- 3 In a fume cupboard, warm but do not boil the acid to help the paper clip to dissolve. The nitric(V) acid oxidises the manganese to $\text{Mn}^{2+}(\text{aq})$ ions.
- 4 Add about 10 cm³ of phosphoric(V) acid to the beaker, followed by about 10 cm³ of potassium iodate(VII) solution. Boil the solution carefully for 10 minutes. Allow the mixture to cool. [The phosphoric(V) acid prevents the precipitation of insoluble iron(III) salts.]
- 5 When the solution is cool pour it into a 100 cm³ volumetric flask using a small funnel. It is important not to lose any of the solution. Rinse the remaining solution from the beaker and funnel into the flask with distilled water and add further distilled water to bring the solution in the flask exactly up to the mark.
- 6 Stopper the flask and shake it to ensure that the solution is uniform. All the manganese that was in the 0.25 g of paper clip is now in the purple solution as the manganate(VII) ion, $\text{MnO}_4^-(\text{aq})$.

(Reference: DENBY, Derek, OTTER, Chris, STEPHENSON, Kay (eds) (2009), *Salters Advanced Chemistry Support Pack*, Heinemann, 71)

I was given a solution of MnO_4^- ions, labelled '2% Mn', which had the same concentration as one that would be produced in my experiment from 0.25 g of steel containing 2% manganese by mass. The solution was actually made up by dissolving potassium manganate(VII) crystals in distilled water in a 1 dm³ volumetric flask. Starting with this solution, I produced the calibration curve (**Fig. 1**) using a colorimeter. I would have preferred to start with a '0.5% Mn' solution, but my teacher said that it wouldn't be as accurate as the '2% Mn' solution.




Fig. 1 % Mn against absorbance

The solution that I made from the paper clips gave an absorbance reading on the colorimeter of 0.64.

Comments on my experiments

I thought that my experiments went well. During the titration I had a bit of a problem with the pipette filler and I think I spilled some of the paper clip solution onto the bench as I was transferring it to the conical flask in one of my titrations.

One extra thing I could have done was to check the concentration of the MnO_4^- ions in the solution that I made from the paper clips by titrating it against a standard solution of Fe^{2+} ions.

END OF PRACTICAL INSERT

Copyright Information:

OCR is committed to seeking permission to reproduce all third-party content that it uses in the assessment materials. OCR has attempted to identify and contact all copyright holders whose work is used in this paper. To avoid the issue of disclosure of answer-related information to candidates, all copyright acknowledgements are reproduced in the OCR Copyright Acknowledgements booklet. This is produced for each series of examinations and is freely available to download from our public website (www.ocr.org.uk) after the live examination series.

If OCR has unwittingly failed to correctly acknowledge or clear any third-party content in this assessment material, OCR will be happy to correct its mistake at the earliest possible opportunity.

For queries or further information please contact the Copyright Team, First Floor, 9 Hills Road, Cambridge CB2 1GE.

OCR is part of the Cambridge Assessment Group; Cambridge Assessment is the brand name of University of Cambridge Local Examinations Syndicate (UCLES), which is itself a department of the University of Cambridge.